

Our Protein Arrays

Our technique: *E. coli* clones are spotted onto PVDF membrane (22 x 22 cm / approx. 8.7 x 8.7 inch), these express recombinant human proteins with N-terminal RGS-His₆ tag.

	hEXselect engine Protein Array	UniPEX engine Protein Array
Product No.	1003	1008
Source	Human fetal brain cDNA	Human fetal brain, T-cells, lung and colon cDNA
Array Content	>57,000 spots in total, >19,000 distinct human antigens, > 4,000 different genes	>36,000 spots in total, >14,000 distinct human antigens, >5,000 different genes

